


# Application Note

## B-01

### Washing Ability Evaluation of Kitchen Detergent



#### Materials

1. Commercial kitchen detergents
2. Protein (Casein: Protein originated from milk.)



#### Results

1. Kitchen detergents could wash the protein dirt immobilized on sensor surface.
2. Washing ability of each detergent was evaluated.



#### Applications

1. Washing evaluation of the effectiveness of various detergents on protein removed. ex.) Kitchen, Contact lens, Bath room, Toilet, Living room and Vehicles.
2. Develop optimum conditions for washes or minimum volume of detergent to be used for washing.
3. Evaluate effectiveness of new detergents.